



visor2™

Precise targeting for TMS therapy

- visor2™ offers a suite of tools that allow TMS providers to precisely target treatment of Major-Depressive Disorder (MDD), Obsessive Compulsive Disorder (OCD), and other” conditions
- Intuitive step-by-step workflows walk through the entire navigated TMS therapy session
- Fast motor threshold determination and precise delivery of TMS pulses
- Simple and precise MRI-less navigation when no individual MRI is available
- Individualized TMS with MRI-assisted neuronavigation and motor threshold determination
- Import of functional and structural data for stimulation target definition (e.g. fMRI)
- Compatible with most TMS coils from Deymed, Neurosoft, Magstim, MagVenture, MAG&More (Apollo), Neuronetics (Neurostar) and Yingchi for cost-effective integration into existing TMS setup

Research over the past several years indicates that with new treatment protocols and targeting parameters the effectiveness of TMS can be dramatically increased (Anish Modak & Paul B. Fitzgerald, 2021; Cole et al., 2020). The visor2™ TMS neuronavigation system positions you to deliver more accurate, consistent, and personalized TMS therapy today, and positions your practice to take advantage of future advancements in TMS targeting.



Personalising transcranial magnetic stimulation for depression using neuroimaging: A systematic review



Stanford Accelerated Intelligent Neuromodulation Therapy for Treatment-Resistant Depression

*visor2 system is CE marked as a medical device in the EU, according to MDD 93/42/EEC, class IIA and has FDA clearance under 510(k) in the USA. Medical Device License (MDL) issued by Health Canada. TGA certified for medical use in Australia.


Manufactured by eemagine GmbH, Berlin, Germany, ISO 13485 certified. ANT Neuro and eemagine are part of the neuromotion group. For the regulatory status of visor2 outside of Australia, Canada, EU and USA, please contact your local distributor or ANT representation. Speech mapping with visor2 is intended for research and educational use only.

The visor2 system is not available for sale in all countries, please refer to the associated visor2 specification documents or contact you local representative for details. <https://www.ant-neuro.com/contact-us>

ANT North America, Inc., Philadelphia, USA,
tel: +1 (608) 204 0878, fax: +1 (608) 204 0874,
internet: www.ant-neuro.com, e-mail: infous@ant-neuro.com

Information in this document is subject to change.

www.ant-neuro.com/products/visor2



High-precision tracking camera

AiO PC with touch-based usability

Mobile medical cart with small footprint



ant neuro
north america